

**KN95 Ohrschlaufen No-Powered Luftreinigungs-
Atemschutzgerät Einweg-Gesichtsschutz**

No-Powered air-purifying respirator

KN95 ear-loops disposable facepiece (dust mask)



— • Specifications • —

Material : non-woven , melt-blown fabric

Color : light blue , White , Black

Size : 108*155mm

Packages : Individual Wrapped Pack , 1 box contains 100pcs
1 carton contains 10 boxes



— • Descriptions • —

Protection guaranteed with over 99% bacterial filtration


Inner Layer

Spun-viscose non-woven, to absorb the hot breath from the body, making the skin dry and comfortable




Mid-Layer

Melt-Blown Fabric, to filter non-oily particles in the air



Outer layer

Spun-viscose non-woven, to effectively block droplet transmission



High elastic ear-loops

with no pressure to the ears

breathable and comfortable



— • Fitting instructions • —


Step1:

Cup the respirator in your hand with the nasal splint up and away from you




Step2:

Position the respirator under your chin and pull the bottom strap around the neck below your ears



Step3:

Pull the top strap over your head resting it high at the top back of your head



Step4:

Place your fingertips of both hands at the middle of the nasal splint. Mould the nasal splint to the shape of your nose. Mould with a single hand may lead to bad seal of the respirator.



— • Take off instructions • —

Step1:

Slowly pull the straps away with your hand. Do not touch the contaminated surface of the mask



Step2:

Fold the mask in half



Step3:

Put the mask in the airtight bag and throw it into the trash bin





Warning !!!

The dust-proof effect will not be achieved by the wrong methods

