**KN95 ear-loops No-Powered air-purifying respirator  
disposable facepiece**

**No-Powered air-purifying respirator**

KN95 ear-loops disposable facepiece (dust mask)

## — • Specifications • —

Material : non-woven , melt-blown fabric

Color : light blue , White , Black

Size : 108\*155mm

Packages : Individual Wrapped Pack , 1 box contains 100pcs  
1 carton contains 10 boxes

# — • Descriptions • —

Protection guaranteed with over 99% bacterial filtration


## Inner Layer

Spun-viscose non-woven, to absorb the hot breath from the body, making the skin dry and comfortable


## Mid-Layer

Melt-Blown Fabric, to filter non-oily particles in the air

## Outer layer

Spun-viscose non-woven, to effectively block droplet transmission

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# — • Fitting instructions • —


## Step1:

Cup the respirator in your hand with the nasal splint up and away from you

## Step2:

Position the respirator under your chin and pull the bottom strap around the neck below your ears

## Step3:

Pull the top strap over your head resting it high at the top back of your head

## Step4:

Place your fingertips of both hands at the middle of the nasal splint. Mould the nasal splint to the shape of your nose. Mould with a single hand may lead to bad seal of the respirator.

# — • Take off instructions • —


## Step1:

Slowly pull the straps away with your hand. Do not touch the contaminated surface of the mask

## Step2:

Fold the mask in half

## Step3:

Put the mask in the airtight bag and throw it into the trash bin

# Warning !!!

The dust-proof effect will not be achieved by the wrong methods