12oz colors mixed glass borosilicate cup

4. Meet FDA test & CA 65 test

1. Different designs and sizes for choice

4. We have exclusive workers for quality control

| Product Details | |
|------------------|---|
| Item No. | SGDS18032703 |
| Material | High borosilicate glass |
| Process | Mouth Blowing Glass |
| Samples time | 1. 5 days if at exist shaped and size of glass |
| | 2. 15 days if need new shape or size of glass |
| Packing | Normal packing,6 pcs into inner box, 72 pcs per carton |
| Product Capacity | 500,000~1,000,000 pcs per month |
| Delivery time | Within 35 days after the sample and order confirmed |
| Payment terms | 30% deposit by T/T in advance and the balance against the copy of B/L |
| Shipment | By sea, by air, by Express and your shipping agent is acceptable |
| | 1. Hand made borosilicate drinking glass _with high quality |
| Product Features | 2. Suitable for used in hotel, home, etc. |
| | 3. It is eco-friendly |

2 .Any color painted, frost, electroplate, pattern laser carver processing3. Special package like shrink wrap, color gift box, white gift box etc.

5. We have professional workshop and warehouse to ensure the delivery time



For your choices











Borosilicate glass



Office & Sample Room



















Method of application:

- 1. Using it under guide of the adult
- 2. Washing it with clean or boiling water before usage
- 3. No touching the rim of glass cup, try to take the bottom or the handle of it

Or here may help you know more about us: FAQ